

1. Product and company identification

Product identifier

Trade name: 636K6 - Plastillin Clay

Relevant identified uses of the substance or mixture and uses advised against

General use: For orthopedic procedures: mass for modeling

Details of the supplier of the safety data sheet

Company name: Otto Bock Health Care

Street/POB-No.: 3820 W. Great Lakes Drive

Postal Code, city: Salt Lake City, UT 84120
USA

WWW: www.ottobockus.com

Telephone: +1 (801) 956-2400

Telefax: +1 (801) 956-2401

Dept. responsible for information:

Quality Department,
Telephone: +1 (801) 954-2304 (7 AM – 3 PM, Mountain Time),
Email: USRegulatory@ottobock.com

Additional information: Corporate headquarters:
Ottobock SE & Co. KGaA
Max-Näder-Straße 15
Duderstadt
Germany

Emergency phone number

CHEMTREC, Telephone: +1 (800) 424-9300

2. Hazards identification

Emergency overview

Appearance: Form: solid

Color: varying colors

Odor: like fatty acid

Classification: This material is classified as not hazardous.

Regulatory status

This material is not considered hazardous by the U.S. OSHA Hazard Communication Standard (29 CFR 1910.1200) and SIMDUT in Canada.

Hazards not otherwise classified

see section 11: Toxicological information

3. Composition / Information on ingredients

Chemical characterization: A mixture of: hydrocarbon waxes, oils, fillers, paraffins, pigment

4. First aid measures

Following skin contact: After contact with molten product, cool skin area rapidly with cold water. Seek medical treatment in case of troubles.

After eye contact: In case of contact with eyes rinse with plenty of water carefully.

Most important symptoms/effects, acute and delayed

After eye contact: mild irritant

Information to physician

Treat symptomatically.

5. Fire fighting measures

Flash point/flash point range:

No data available

Auto-ignition temperature: No data available

Suitable extinguishing media:

Water spray jet, foam, dry chemical powder, carbon dioxide.

Extinguishing media which must not be used for safety reasons:

Full water jet

Specific hazards arising from the chemical

In case of fire may be liberated: carbon monoxide and carbon dioxide

Protective equipment and precautions for firefighters:

Wear self-contained breathing apparatus.

6. Accidental release measures

Personal precautions: Handle in accordance with good industrial hygiene and safety practice.

Environmental precautions:

Discharge into the environment must be avoided.

Methods for clean-up:

Take up mechanically, placing in appropriate containers for disposal.

7. Handling and storage

Handling

Advices on safe handling: Handle in accordance with good industrial hygiene and safety practice.

Specific use(s)

For orthopedic procedures: mass for modeling

Storage

Requirements for storerooms and containers:

Keep container tightly closed and dry. Keep in a cool place.

8. Exposure controls / personal protection

Engineering controls

Provide adequate ventilation, and local exhaust as needed.

See also information in chapter 7, section storage.

Personal protection equipment (PPE)

Eye/face protection	Wear safety goggles when handling hot molten mass. According to OSHA Standard - 29 CFR: 1910.133 or ANSI Z87.1-2003.
Skin protection	If necessary: protective gloves according to OSHA Standard - 29 CFR: 1910.138 Glove material: nitrile rubber-breakthrough time: 480 min. Observe glove manufacturer's instructions concerning penetrability and breakthrough time.
General hygiene considerations:	Protect from excessive heat. Wash hands before breaks and after work.

9. Physical and chemical properties

Information on basic physical and chemical properties

Appearance:	Form: solid Color: varying colors
Odor:	like fatty acid
Odor threshold:	No data available
pH value:	No data available
Melting point/freezing point:	No data available
Initial boiling point and boiling range:	No data available
Flash point/flash point range:	No data available
Evaporation rate:	No data available
Flammability:	No data available
Explosion limits:	No data available
Vapor pressure:	No data available
Vapor density:	No data available
Density:	1.60 g/cm ³
Water solubility:	insoluble
Partition coefficient: n-octanol/water:	No data available
Auto-ignition temperature:	No data available
Thermal decomposition:	approx. 536 °F
Ignition temperature:	approx. 662 °F

10. Stability and reactivity

Reactivity:	refer to section 10.3
Chemical stability:	Product is stable under normal storage conditions.
Possibility of hazardous reactions	No dangerous reactions are known.
Conditions to avoid:	Protect from excessive heat.
Incompatible materials:	No data available
Hazardous decomposition products:	In case of fire may be liberated: carbon monoxide and carbon dioxide
Thermal decomposition:	approx. 536 °F

11. Toxicological information

Toxicological tests

Toxicological effects: Acute toxicity (oral): Lack of data.
 Acute toxicity (dermal): Lack of data.
 Acute toxicity (inhalative): Lack of data.
 Skin corrosion/irritation: Lack of data.
 Serious eye damage/irritation: Lack of data.
 Sensitisation to the respiratory tract: Lack of data.
 Skin sensitisation: Lack of data.
 Germ cell mutagenicity/Genotoxicity: Lack of data.
 Carcinogenicity: Lack of data.
 Reproductive toxicity: Lack of data.
 Effects on or via lactation: Lack of data.
 Specific target organ toxicity (single exposure): Lack of data.
 Specific target organ toxicity (repeated exposure): Lack of data.
 Aspiration hazard: Lack of data.

Other information: In case of heating: risk of burns.

Symptoms

After eye contact: mild irritant

12. Ecological information

Ecotoxicity

Further details: No data available

Mobility in soil

No data available

Persistence and degradability

Further details: No data available

Additional ecological information

General information: Discharge into the environment must be avoided.

13. Disposal considerations

Product

Recommendation: Incinerate according to applicable local, state and federal regulations.

Contaminated packaging

Recommendation: Dispose of waste according to applicable legislation.
 Non-contaminated packages may be recycled.

14. Transport information

USA: Department of Transportation (DOT)

Proper shipping name: Not restricted

Sea transport (IMDG)

Proper shipping name: Not restricted

Marine pollutant: no

Air transport (IATA)

Proper shipping name: Not restricted

Further information

No dangerous good in sense of these transport regulations.

15. Regulatory information

National regulations - Great Britain

Hazchem-Code: -

16. Other information

Hazard rating systems:



NFPA Hazard Rating:

Health: 0 (Minimal)

Fire: 1 (Slight)

Reactivity: 0 (Minimal)

HMIS Version III Rating:

Health: 0 (Minimal)

Flammability: 1 (Slight)

Physical Hazard: 0 (Minimal)

Personal Protection: X = Consult your supervisor

HEALTH	0
FLAMMABILITY	1
PHYSICAL HAZARD	0
	X

Reason of change: Changes in section 1.3: Corporate headquarters

Date of first version: 10/15/2008

Department issuing data sheet

Contact person: see section 1: Dept. responsible for information

The information in this data sheet has been established to our best knowledge and was up-to-date at time of revision. It does not represent a guarantee for the properties of the product described in terms of the legal warranty regulations.