


636W28 special adhesive and 756K30 dispensing gun

Dispensing and gluing made easy

Based on epoxy resin, 636W28 adhesive is suitable for applications including gluing in foot adapters as well as bonding wood and metal, unilateral bar systems and E-MAG knee joints.

For convenient processing, we offer the 756K30 dispensing gun for the 50 ml container that makes it easy to squeeze the special adhesive out of cartridges A and B.

Thanks to the enclosed static mixers, the adhesive can be used multiple times. This way, the adhesive doesn't dry too quickly and is still ready to use even over a longer period of time.



Benefits at a glance:

- 636W28 special adhesive is suitable for gluing in foot adapters as well as bonding wood and metal, unilateral bar systems and E-MAG knee joints
- 756K30 dispensing gun makes it easier to dispense the special adhesive with the 50 ml container
- Can be used multiple times as the adhesive does not dry out immediately

636W28 Special adhesive

Article number 636W28

Net contents 12.5 ml adhesive and 12.5 ml hardener

Package contents One cartridge, one plunger and five static mixers

Article number 636W28=0.050

Net contents 25 ml adhesive and 25 ml hardener

Package contents One cartridge and five static mixers

Article number 756K30

Package contents One dispensing gun