

Ottobock plaster modelling aid

Your trusted helper for plaster working



Quality for life

Shaping plaster models quickly and easily

With the plaster modelling aid, fabrication of plaster models is easier than ever before. Continuous clamping and reclamping of the plaster model is now a thing of the past. You save valuable time and can modify your models more efficiently.

Here's how it works:

- 1.** Clamp the plaster modelling aid into the vice.
- 2.** Secure the plaster model in the clamping device.
- 3.** Begin modifying the plaster model. During modelling you can rotate and reorient the workpiece at any time – easily and with just one hand.

As you can see, with this innovative plaster modelling aid you can modify plaster models quickly and efficiently. The advantages to you are manifold: you can optimise your fabrication process, save time and make your daily work easier.

**Product characteristics
plaster modelling aid
article no. 704G7**

- Sturdy design
- Easy handling
- Weight: 1.7 kg
- Material: stainless steel
- Dimensions (L x W x H):
119 mm x 55 mm x 130 mm
- Suitable for tube diameters of
16 mm, 21 mm, 26 mm and 36 mm

