

# Carbon Ankle Seven

## iFab Order Form



Account Information		Ship To Information	
Date	<input type="text"/>	Name	<input type="text"/>
Account Number	<input type="text"/>	Address	<input type="text"/>
Bill to	<input type="text"/>	City/State/Zip	<input type="text"/>
Phone Number	<input type="text"/>	Phone Number	<input type="text"/>
Email	<input type="text"/>	Email	<input type="text"/>
Buyer	<input type="text"/>	Desired Delivery Date	<input type="text"/>
PO Number	<input type="text"/>		

Please mail your completed Ottobock order form and a negative impression of the patient's limb to the address below. An Ottobock Fabrication Coordinator will contact you.

**Shipping Options:**

UPS Next Day      UPS Ground      UPS 2-Day      Other

Once this form is complete, you may submit via email, or print and include in the shipment with the patient's negative impression.

**Order      Quote Only**

**For clinical questions, call 800 328 4058.**

**Notes:**



**Patient Information:**

Left

Patient Name \_\_\_\_\_ Patient Weight \_\_\_\_\_ lbs. Right Pathology \_\_\_\_\_

**Select Type of AFO**

Removable custom anterior tibial shell  
Split-front floor reaction

**Material Selection**

Copoly  
Specify Friddle's Pattern\* (please specify below)

\*Additional charges may apply for Friddle's patterns.

**Padding Material**

Microcell Puff  
Specify color:  White  Black

Bocklite (Only available in beige)

**Material Color**

White  
Black

**Plastic Color**

White  
Black

**CARBON ANKLE SEVEN WEIGHT & ACTIVITY SELECTION GUIDE**

Body Weight	Activity Level		Spring Width
	Normal	High	
220lbs	17CF1=L/R1	17CF1=L/R1	30mm
198lbs	17CF1=L/R2	17CF1=L/R1	30mm
176lbs	17CF1=L/R3	17CF1=L/R2	30mm
154lbs	17CF1=L/R4	17CF1=L/R3	30mm
132lbs	17CF1=L/R5	17CF1=L/R4	30mm
110lbs	17CF1=L/R6	17CF1=L/R5	25/30mm
88lbs	17CF1=L/R8	17CF1=L/R7	25mm
66lbs	17CF1=L/R9	17CF1=L/R8	25mm
44lbs	17CF1=L/R11	17CF1=L/R10	22mm
22lbs	17CF1=L/R12	17CF1=L/R11	22mm

**Correction to Cast**

Correct to match Carbon Ankle Seven strut  
Leave as casted  
to 90° Valgus/varus correction  
to \_\_\_\_\_ ° plantarflexion Tone reducing feature  
to \_\_\_\_\_ ° dorsiflexion

**Other Options**

All Carbon Ankle Sevens come with a full liner unless otherwise indicated.

Malleolar pads Heel post \_\_\_\_\_ mm

Medial arch support pad

Standard Carbon Ankle Seven AFO supplied with one 2-inch hook-and-loop strap and chafe fastened on the lateral side.

Check here if you would like to use the heel height of the Carbon Ankle Seven strut, which is the recommended heel height.

**Metric Measurements Only**

