

# Genium X4 Cover Solutions.

Explore the comprehensive range of options and configurations for our next-gen microprocessor knee (MPK).



## Select the right functionality

See for yourself



- Choose from a variety of covers, optimized for different activities and environments.
- Pick the functional, protective, or cosmetic option that best suits specific needs, goals, and preferences.

## Choose the optimal components

- Shield specific areas, or select a cover that provides top-to-bottom protection – as well as a natural look, if desired.
- Enhance comfort and confidence without adding significant weight.

## Create the perfect look

See for yourself



- Many styling options can be personalized with decorative elements, such as spray paint and stickers.
- Create a combination that lets individual style and personality shine through.

## Genium X4 Cover Solutions.

Select and customize the cover that provides exactly the right function, protection, and style.



### Protective Cover (Short)

**Article No:** 4P100=7  
**Care & Handling:** Easily removable, dishwasher-safe, personalizable

**Details:** Made from durable thermoplastic polyurethane (TPU) material, this sleek and sporty cover solution protects the *Genium X4* from regular wear and tear as well as tough, rugged conditions. Plus, it can be decorated with spray paint, stickers, and more based on personal style and aesthetic preferences.



### Functional Knee Part

**Article No:** 4P112=1  
**Care & Handling:** Easily removable, dishwasher safe

**Details:** The functional knee part (FKP) acts like a kneecap, giving the *Genium X4* a more natural appearance without restricting movement. It can be worn on it's own, as shown here, or combined with the tube cover.



### Protective Cover (Long)

**Article No:** 4P110=7  
**Care & Handling (All):** Easily removable  
**(Textile Cover):** Machine washable, personalizable  
**(Tube Cover):** Dishwasher safe

**Details:** The long protective cover consists of the tube cover (see more info below) wrapped in a durable, personalizable TPU and polyurethane textile cover. It keeps the *Genium X4* knee joint and tube adapter shining like new amidst all the bumps, knocks, scrapes, dust, dirt, and more that keep life exciting.



### Tube Cover\*

**Article No:** \*Available only as part of protective cover (long)  
**Care & Handling:** Easily removable, dishwasher safe

**Details:** The tube cover includes a tube sleeve and tube clamp. It not only protects the tube adapter from dirt and scratches, but it also provides a smooth transition from the knee down to the foot. Although it's only available as part of the long protective cover solution (shown above), it can be worn on its own or combined with the functional knee part (4P112=1).



### Functional Cosmesis

**Article No (knee part + shank):** 3F2=0  
**Article No (thigh kit):** 3D13=1  
**Article No (stocking):** 99B122 (size: S, L, color: beige (-4), light brown (-15))  
**Care & Handling (All):** Easily removable  
**(Stocking):** Water repellent, machine washable

**Details:** This cover solution provides the ultimate natural look without restricting movement. It combines a knee part, a foam shank that can be formed to match the contralateral leg, a stocking available in light and dark colors, and an optional thigh kit with magnetic closure to smooth the transition from thigh to cosmesis.