

## FS3 Defender

### Product features

- For children and smaller users
- Ensures a smooth rollover and exceptional ground compliance in a flexible design
- Designed to withstand the particularly high activity level of children



### Technical data

K-Level	3, 4
Max. body weight	100 lbs / 45 kg
Sizes	19-21 cm
Heel height	3/8" / 10 mm
Weight with footshell*	7.94 oz / 225 g
Build height*	3" / 76 mm

\* Reference size is 20 cm.

## FS3 Defender

### Order form

Customer		Shipping address (if different from customer address)	
Customer no.	<input type="text"/>	Customer no.	<input type="text"/>
Company	<input type="text"/>	Company	<input type="text"/>
Street	<input type="text"/>	Street	<input type="text"/>
Postal code/city	<input type="text"/>	Postal code/city	<input type="text"/>
O&P professional	<input type="text"/>	Customer ref.	<input type="text"/>

#### Stiffness chart

Body weight	Stiffness
20–30 lbs / 9–13 kg	1P
31–40 lbs / 14–18 kg	2P
41–50 lbs / 19–22 kg	3P
51–60 lbs / 23–27 kg	4P
61–70 lbs / 28–31 kg	5P
71–80 lbs / 32–36 kg	6P
81–100 lbs / 37–45 kg	1A (Adult)

#### Scope of delivery

- FS3\* Defender
- KIT-00-12201-00 Set of heel wedges
- SO-NPS-200\* Spectra protective sock black

#### Accessories & spare parts

- SO-NPS-100\* Spectra protective sock white

#### Available footshell variants

- FTC-2M-1\*: Defender footshell

#### Foot module part number

Reference number	Stiffness	Size
FS3 - 00 - 0		RU
	Stiffness	Size (cm)
	1P–6P	19, 20, 21
	1A	

#### Footshell part number

Reference number	Size	Color	Side
FTC - 2M - 1	4		R
	Size (cm)	Color	Side
	19, 20, 21	L Light	L Left
		M Medium	R Right
		D Dark	