




## More than a foot. A foundation

Each of your patients is different. They all move to their own rhythm and speed. They have individual habits and preferences, and their own goals.

That is why we at Ottobock have continued to develop our portfolio over the last 100 years. For us, mobility is more than a level or a unit of measurement – it's an attitude to life. We want to provide the best possible support for users, step by step.

And our development continues: now with Taleo, Taleo Low Profile, Taleo Vertical Shock, Taleo Harmony, and Taleo Side Flex, our expanded family of feet provides you with more options to choose the right foot for your patient.



## ottobock. The Ottobock Foot Portfolio

Quality for life

Information for practitioners

Ottobock US  
P 800 328 4058  
F 800 962 2549  
professionals.ottobockus.com

Ottobock Canada  
P 800 665 3327  
F 800 463 3659  
professionals.ottobock.ca

©2021 Otto Bock HealthCare LP · 19616 · 07/21



## Terion K2

Confidence in every step.



Designed for less active individuals who generally navigate indoor environments and utilize a walking aid for outdoor ambulation and place a high value on dependable support from their prosthetic foot.

### Benefits of 1C11 Terion K2:

- Comfortable heel strike and fast full ground contact for improved stability
- Easy rollover at slow walking speeds
- Increased confidence due to stability while standing and walking

K1, K2 Up to 385 lbs / 175 kg



## Trias

Secure as expected.



Designed for moderately active individuals who navigate indoor and familiar outdoor environments and place a high value on consistent stability when walking.

### Benefits of 1C30-1 Trias:

- Soft rollover for easy movement at slow to moderate walking speeds
- The right amount of energy for moderate activities
- Supports controlled movements across stable terrain

K3 Up to 275 lbs / 125 kg Waterproof: fresh, salt and chlorinated



## Taleo

Ready for everyday life.



Designed for active individuals who navigate varied indoor and outdoor environments and place a high value on effortless walking and the ability to go wherever life takes them.

### Benefits of 1C5\* Taleo:

- Smooth rollover for effortless walking
- Energy efficient walking for a variety of activities
- Supports dynamic movement across varying terrain conditions

K3 Up to 330 lbs / 150 kg Waterproof: fresh, salt and chlorinated<sup>1</sup>

\* Family members: 1C50 Taleo<sup>1</sup>, 1C51 Taleo Vertical Shock, 1C52 Taleo Harmony, 1C53 Taleo Low Profile<sup>1</sup>, 1C58 Taleo Side Flex<sup>1</sup>.



## Triton

Your will. Your way.



Designed for highly active individuals who navigate varied indoor and outdoor environments and place a high value on uncompromised response and control even when performing high impact activities.

### Benefits of 1C6\* Triton:

- Dynamic rollover supports vigorous walking
- High energy return for every adventure
- Supports agile and high-impact movement across varying terrain conditions

K3, K4 Up to 330 lbs / 150 kg Waterproof: fresh, salt and chlorinated<sup>1</sup>

\* Family members: 1C60 Triton, 1C61 Triton Vertical Shock, 1C62 Triton Harmony, 1C63 Triton Low Profile<sup>1</sup>, 1C64 Triton Heavy Duty<sup>1</sup>, 1C68 Triton side flex<sup>1</sup>.