


# Hand Portfolio Selection Guide




## Thumb



## Wrist



## Wrist & Thumb



## Children



# Hand Orthoses Overview

## Thumb


50P22 Rhizo Arexa



- Formable lateral aluminium splints
- No pressure on the back of the thumb

Side	Size	Circumference (cm)
Bilateral	S	14–18
	M	18–22

4085 Thumboform Long



- Patient-specific fit
- Formable thermoplastic plastic insert

Side	Size	Circumference (cm)
L/R	XS	13–14
	S	14–16
L/R	M	16–18
	L	18–20
L/R	XL	20–22

## Wrist

50P70 Manu Carezza



- Sporty design for activities
- Handy closure strap

Size	Circumference (cm)
Universal	14–22

50P13 Manu Sensa



- Three-dimensional flat-knitted fabric
- Anatomically pre-shaped volar plastic splint

Side	Size	Circumference (cm)
L/R	XS	14–15
	S	15–16
L/R	M	16–17
	L	17–18
L/R	XL	18–19
	XXL	19–20

• Colors: 1= grey; 0= beige

4142 Manu 3D



- Anatomically pre-shaped volar plastic splint
- Handy closure straps, e.g. for swollen wrists

Side	Size	Circumference (cm)
L/R	XS	14–15
	S	15–16
L/R	M	16–17
	L	17–18.5
L/R	XL	18.5–20
	XXL	20–22

• Colors: 1= grey; 7= black

4143 Manu 3D Basic




- Anatomically pre-shaped volar plastic splint
- Adjustable thumb loop

Side	Size	Circumference (cm)
L/R	XS	14–15
	S	15–16
L/R	M	16–17
	L	17–18.5
L/R	XL	18.5–20
	XXL	20–22

• Colors: 1= grey; 7= black

50P20 Manu Arexa



- Anatomically pre-shaped volar, dorsal, or ulnar aluminium splint
- Pull-2-Go closure for easy, one-handed use

Side	Size	Circumference (cm)
Bilateral	S	14–16
	M	16–19
	L	19–22

## Wrist


4088 Wristoform



- Patient-specific fit
- Formable thermoplastic volar splint and plastic dorsal insert

Side	Size	Circumference (cm)
L/R	XS	13–14
	S	14–16
L/R	M	16–18
	L	18–20
L/R	XL	20–22

50P10 Manu Immobil



- Individually adjustable hand position
- Washable materials for post-operative use

Size	Hand Length (cm)
Universal	< 18
	> 18

50P11 Manu Immobil Long



- Individually adjustable hand position
- Additional individually adjustable finger support

Size	Hand Length (cm)
M	< 18
L	> 18

## Wrist & Thumb

4146 Manu 3D Pollex



- Anatomically pre-shaped volar plastic splint and coil spring for the thumb
- Handy closure straps, e.g. for swollen wrists

Side	Size	Circumference (cm)
L/R	XS	14–15
	S	15–16
L/R	M	16–17
	L	17–18.5
L/R	XL	18.5–20
	XXL	20–22

• Colors: 1= grey; 7= black

50P21 Manu Arexa Pollex



- Anatomically pre-shaped volar and dorsal aluminium splint and coil spring for the thumb
- Pull-2-Go closure for easy, one-handed use

Side	Size	Circumference (cm)
Bilateral	S	14–16
	M	16–19
	L	19–22

## Children

50P50 Manu Junior



- Easy to use
- Anatomically pre-shaped volar aluminium splint

Side	Size	Circumference (cm)
L/R	S	9–11
	M	11–13
L/R	L	13–15

50P51 Manu Junior Pollex



- Easy to use
- Anatomically pre-shaped volar aluminium splint and additional support for the thumb

Side	Size	Circumference (cm)
L/R	S	9–11
	M	11–13
L/R	L	13–15