ottobock.

LP2-W2 Freestyle Swim

Product features

- With the Freestyle Swim, users can walk comfortably to and from the swim area
- For swimming with flippers, only the button has to be pressed to achieve a high swimming speed, which reduces energy expenditure
- The ankle joint can be fixed in 0° or with 70° plantar flexion
- The prosthetic foot supports the balancing of the body, enabling efficient swimming



Technical data

K-Level	3, 4
Max. body weight	255 lbs / 116 kg
Sizes	22-31 cm
Heel height	3/8" / 10 mm
Weight with footshell*	30.34 oz / 860 g
Build height*	4 ½" / 103 mm
••••••	• • • • • • • • • • • • • • • • • • • •

^{*} Reference size is 26 cm.

ottobock.

LP2-W2 Freestyle Swim

Order form

	Customer		Shipping address (if different from customer address)
Customer no.		Customer no.	
Company		Company	
Street		Street	
Postal code/city		Postal code/city	
0&	P professional		Customer ref.

Stiffness chart

Body weight		Activity level	
	Low	Moderate	High
100-115 lbs / 44-52kg	1	1	2
116-130 lbs / 53-59 kg	1	2	3
131-150 lbs / 60-68 kg	2	3	4
151-170 lbs / 69-77 kg	3	4	5
171–195 lbs / 78–88 kg	4	5	6
196-220 lbs / 89-100 kg	5	6	_
221-255 lbs / 101-116 kg	6	_	_

Foot module part number

Reference number					Stiffness		Size			
LP2 -		W2	-	0	- '					RU
					L		₹	*		•••••
	Stiffness			Size (cm)						
						1-6		22, 23,, 31		

Footshell part number

Reference number					Size	Color					Side
FTC	-	3C	-	1				4	-	R	
					\	• • • • • • • • • • • • • • • • • • • •	L	▼			Ц,
					Size (cm)		Color				Side
					22, 23,, 31		L Light				L Left
					•	••••••••••••••••••	M Medi	um			R Right
							D Dark				••••••

Scope of delivery

- LP2-W2* Freestyle Swim
- KIT-00-11* Set of heel wedges
- S0-NPS-200* Spectra protective sock black

Accessories & spare parts

- S0-NPS-100* Spectra protective sock white
- ACC-00-10300-00 Footshell removal tool
- KIT-00-20000-00 Protective ankle cover

Available footshell variants

• FTC-3C-1*: Footshell Freestyle Swim (without cap)