





# MyoBock Hand Option Chart

	Digital Twin (*=7)	DMC Plus (*=6)	VariPlus Speed (*=9)	SensorHand Speed (*=8)
<b>Proportional Control</b>	✗	✓	✓	✓
<b>Opening Speed</b>	110 mm/sec	15-130 mm/sec	15-300 mm/sec	15-300 mm/sec
<b>Grip Strength</b>	90 N (20 lb)	0-90 N (0-20 lb)	0-100 N (0-22 lb)	0-100 N (0-22 lb)
<b>Coaxial Receiver Color</b>	 White	 Brown	 Green	 Orange
<b>Programming</b>	Function Plug	Function Plug	MyoSelect	MyoSelect or Function Plug
<b>Electrode Control</b>	Single or Dual Site	Dual Site Only	Single or Dual Site	Single or Dual Site
<b>Myo Plus TR Pattern Recognition Compatible</b>	✗	✗	✓	✓
<b>Quick Disconnect Wrist (8E38*)</b>	Sizes: 7 ¼, 7 ¾, 8 ¼	Size 7 Only	Sizes: 7 ¼, 7 ¾, 8 ¼	Sizes: 7 ¼, 7 ¾, 8 ¼
<b>Wrist Disarticulation (8E39*)</b>	Sizes: 7 ¼, 7 ¾, 8 ¼	Size 7 only	Sizes: 7 ¼, 7 ¾, 8 ¼	Sizes: 7 ¼, 7 ¾, 8 ¼
<b>Transcarpal Wrist (8E44*)</b>	✗	Sizes: 7 ¼, 7 ¾, 8 ¼	✗	✗
<b>Threaded Stud (8E41*)</b>	✗	✗	Sizes: 7 ¼, 7 ¾, 8 ¼ Optional Flexion Wrist (USA only through MASC)	✗