

VS2 Promenade

Product features

- Near effortless rollover optimized for more comfort
- Partial split-toe offers controlled inversion/eversion, providing users enhanced support on uneven surfaces
- Full-length sole plate facilitates a smooth transition from heel strike all the way through toe off, providing a more fluid gait and increased comfort throughout the day



Technical data

K-Level	3
Max. body weight	325 lbs / 147 kg
Sizes	22-31 cm
Heel height	3/8" / 10 mm
Weight with footshell*	20.99 oz / 595 g
Build height*	4 7/8" / 124 mm

* Reference size is 26 cm.

VS2 Promenade

Order form

Customer		Shipping address (if different from customer address)	
Customer no.	<input type="text"/>	Customer no.	<input type="text"/>
Company	<input type="text"/>	Company	<input type="text"/>
Street	<input type="text"/>	Street	<input type="text"/>
Postal code/city	<input type="text"/>	Postal code/city	<input type="text"/>
O&P professional	<input type="text"/>	Customer ref.	<input type="text"/>

Stiffness chart

Body weight	Stiffness
100–126 lbs / 44–57 kg	1
127–160 lbs / 58–72 kg	2
161–203 lbs / 73–92 kg	3
204–258 lbs / 93–117 kg	4
259–325 lbs / 118–147 kg	5

Scope of delivery

- VS2* Promenade
- KIT-00-113* Set of heel wedges
- SO-NPS-200* Spectra protective sock black

Accessories & spare parts

- SO-NPS-100* Spectra protective sock white
- ACC-00-10300-00 Footshell removal tool

Foot module part number

Reference number	Stiffness	Size
VS2 - 00 - 0	A	- RU
	Stiffness 1–5	Size (cm) 22, 23, ..., 31

Available footshell variants

- FTC-3M-0*: Footshell standard (with cap)
- FTC-3M-1*: Footshell standard (without cap)

Footshell part number

Reference number	Generation	Profile	Cap	Size	Color	Side
FTC -				4		R
	Generation 2–3	Profile M Mid H High	Cap 0 Yes 1 No	Size (cm) 22, 23, ..., 31	Color L Light M Medium D Dark	Side L Left R Right